



Amplicon Sequencing



Amplicon sequencing is frequently used to identify and differentiate microbial species. Short (< 470 bp) hypervariable regions of conserved genes or intergenic regions are amplified by PCR, analyzed using NGS technology, and the resulting sequences are compared against microbial databases.


For bacteria and archaea, the 16S rRNA gene is the most common target for amplicon sequencing. For fungi, three targets are generally used: the 18S rRNA gene, and two internal transcribed spacers (ITS) located between rRNA genes. These regions are usually sufficiently divergent to separate even highly related species, and can sometimes differentiate subspecies.

At Medikonia, we have sequenced over 50,000 microbial samples for our customers. Our standard bioinformatics analyses include OTUs cluster, species annotation, alpha-diversity analysis, beta-diversity analysis, and multi-variate statistical analysis. Applications range from identifying a single species in pure culture to characterizing the microbiota of animals or plants to comparing species diversity and population structure in various environmental sources or geographic regions. Our specialists can advise you on the appropriate analyses for your project.

The Advantages

- **Highly experienced:** We have sequenced over 50,000 samples, resulting in nearly 20 published articles.
- **Outstanding service:** We provide high-quality sequencing, an efficient standard workflow, fast turnaround time, and bioinformatics analyses at a cost-effective price.
- **Effective methodology:** Our method features high amplification efficiency of sample DNA (>95%) and uses PCR free libraries to avoid amplification bias.
- **Comprehensive analysis:** We provide expert bioinformatics analyses using the latest sequence databases and software, generating high-quality, publication-ready data.

Project Workflow



SEQUENCING STRATEGY

- 130-470bp insert DNA library
- HiSeq platform, paired-end 250bp

DATA QUALITY GUARANTEE

- The amount of data for each sample is not less than 30,000 tags, 50,000 tags or 100,000 tags.

SAMPLE REQUIREMENTS

- DNA amount: ≥ 150 ng (for one library preparation*)
- Multiple DNA samples can be mixed together to make one library construction.
- PCR product ≥ 400 ng (for one library preparation)
- DNA concentration: ≥ 20 ng/ μ l
- Sample Volume: >20 μ l
- Purity: OD260/280 = 1.8 - 2.0 without degradation and RNA contamination
- Amplified region should be less than 470 bp

TURNAROUND TIME


- Within 15 working days from verification of sample quality (without data analysis)
- Additional 11 working days for data analysis

RECOMMENDED SEQUENCING DEPTH

- Three strategies: 30,000 tags, 50,000 tags, or 100,000 tags

Target	Region	Fragment Length	Primer	Primer Sequences (5'-3')
Bacterial 16S rDNA	V4	292 bp	515F	GTGCCAGCMGCGCGGTAA
			806R	GGACTACHVGGGTWTCTAAT
	V3-V4	466 bp	341F	CCTAYGGGRBGCASCAG
			806R	GGACTACNNGGTATCTAAT
	V4-V5	393 bp	515F	GTGCCAGCMGCGCGGTAA
			907R	CCGTCAATTCCTTTGAGTTT
Archaeal 16S rDNA	V4	288 bp	U519F	CAGYMGCCRCGGAHAACC
			806R	GGACTACNSGGGTMTCTAAT
	V4	397bp	519F	CAGCCGCCGCGGTAA
			915R	GTGCTCCCCCGCAATTCCT
Eukaryotic 18S rDNA	V4	179 bp	528F	GCGGTAATTCCAGCTCCAA
			706R	AATCCRAGAATTCACCTCT
	V9	131 bp	1380F	CCCTGCCHTTTGTACACAC
			1510R	CCTTCYGCAGGTTACCTAC
Fungal ITS*	ITS1	307 bp	ITS5-1737F	GGAAGTAAAAGTCGTAACAAGG
			ITS2-2043R	GCTGCGTTCTTCATCGATGC
	ITS2	386 bp	ITS3	GCATCGATGAAGAACGCAGC
			ITS4	TCCTCCGCTTATTGATATGC

Analysis Pipeline



LIST OF ANALYSES

- Data quality control
- OTUs cluster and phylogenetic relationship construction
- Species annotation
- Alpha diversity analysis (Observed species, Goods coverage, Chao1, ACE, Shannon, Simpson Index)
- Beta diversity analysis (Unifrac distance heatmap, PCA, PCoA, UPGMA)
- Multi-variate statistical analysis, including NMDS analysis, LEfSe (LDA effect size) analysis, metastats analysis, species T-test analysis, Anodim and MRPP analysis.

Project Example


The following study utilized Medikonia's amplicon sequencing services.

Microbiota modulate tumoral immune surveillance in lung through a $\gamma\delta T17$ immune cell-dependent mechanism

Cancer Research, 74:4030 (2014)

In this study on the effects of commensal bacteria on immune homeostasis, Medikonia technology was used to monitor the changes in commensal bacteria following antibiotic treatment, using stool samples from treated and untreated mice. Analysis of 16S rDNA using the Illumina MiSeq platform showed that multiple taxonomic groups were present, and that their relative abundance was altered by antibiotic treatment.

Microbiota variations after different antibiotic treatments



EXAMPLES OF PUBLICATIONS

Journal	Title
Cancer Research, 74:4030 (2014)	Microbiota modulate tumoral immune surveillance in lung through a $\gamma\delta T17$ immune cell-dependent mechanism.
Scientific Reports, 5:9342 (2015)	The bacterial communities associated with fecal types and body weight of rex rabbits.
Environmental science & technology, 49 (12):7152 (2015)	Community Structure and Soil pH Determine Chemoautotrophic Carbon Dioxide Fixation in Drained Paddy Soils.
Journal of Hazardous Materials, 293:37 (2015)	Conductive iron oxide minerals accelerate syntrophic cooperation in methanogenic benzoate degradation
Journal of hazardous materials, 286:457 (2015)	Exploring bacterial community structure and function associated with atrazine biodegradation in repeatedly treated soils.
Bioresource Technology, 167:124 (2014)	Enrichment of anodic biofilm inoculated with anaerobic or aerobic sludge in single chambered air-cathode microbial fuel cells.